


OMNI2 - SIMULATION AND DEBRIEFING CONTROLLER

PRODUCT NUMBER: 1706



DESCRIPTION:

OMNI 2 makes it easier than ever to drive scenarios, monitor performance, and capture participants' actions for data-rich debriefing sessions.

It is simple to operate with touchscreen controls and an intuitive layout that lets you manage physiological changes while remaining focused on the training. Put simply, it's a frustration-free solution that works to complement your clinical know-how.

OMNI 2 is compatible with a growing selection of over 40 Gaumard patient simulators and skills trainers.

CASES AND PATHOLOGIES:

NURSING SKILLS:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- ▣ Simulate realistic clinical presentations to train patient assessment and care
- ▣ Train patient monitor interpretation and management planning
- ▣ Features skill-specific modules including blood pressure and CPR trainers with real time feedback

OBSTETRIC SKILLS:

- ▣ Automatic fetus delivery controls offer you repeatability and consistency between exercises
- ▣ Rehearse management of shoulder dystocia
- ▣ Virtual CTG monitor support offers practice in FHR/UA interpretation, decision making and documentation

NEONATAL RESUSCITATION:

- ▣ Simulate respiratory compromise to train assessment and management protocols
- ▣ Monitor ventilation performance in real-time to improve technique
- ▣ Built-in resuscitation algorithms help track participant actions and key events for debriefing

ECPR:

- ▣ Real-time CPR feedback: Compression depth and rate, excessive ventilation, no-flow time, CPR cycle
- ▣ Smart Coach: audible tones make it easy to teach the correct rate and ratio every time
- ▣ Detailed report lets you identify and improve the weak areas in a provider's performance

CHARACTERISTICS:

CONTROL:

- ▣ Wireless communication lets you move freely to better observe the action
- ▣ Over 35 programmable vitals including: HR, ECG, RR, BP, SpO2, EtCO2, and more
- ▣ One-touch vitals controls facilitate changes on-the-fly or trending over time

MONITORING:

- ▣ Monitor CPR quality metrics in real time to ensure quality
- ▣ Supports optional virtual patient monitor for training assessment and decision-making skills
- ▣ Monitor ventilation quality and its effectiveness on the patient

DEBRIEFING:

- ▣ Timestamped log records provider actions, vital signs changes, and notes to aid debriefing
- ▣ Includes preprogrammed lists of common provider actions for easy tracking
- ▣ Save and share session log for archiving and debriefing

COMPATIBLE PATIENT SIMULATORS:

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- ☒- Code Blue III family
- ☒- HAL family
- ☒- SUSIE and CHLOE family
- ☒- NOELLE family
- ☒- PREMIE Blue, PEDI Blue family
- ☒- SUSIE SIMON simulators

COMPATIBLE SKILLS TRAINERS:

- ☒- Blood pressure reading skills trainer
- ☒- HAL heart and lung sounds family

INCLUDE:

- ☒- 1 OMNI2 tablet controller
- ☒- 1 Rugged case
- ☒- 1 Power supply + cable
- ☒- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>