

Advanced TLS and PHTLS simulator

PRODUCT NUMBER: 1342



DESCRIPTION:

This Packable Wound Trainer simulates the procedures of **intraosseous insertions, needle decompression and wound packing.**

The realistic wound is located by the hip and can be filled with blood via a bleeding tube. The two wounds include a **rear gunshot entry and a deep frontal exit wound** with a large internal cavity (roughly the size of a golf ball). These include built-in tubes with bleeding capabilities.

SKILLS:

- ☑- Wounds to provide training on wound packing
- ☑- Intraosseous insertions in bilateral humeral head
- ☑- Needle decompression: 2nd intercostal space in the midclavicular line (front)
- ☑- Includes bleeding tubes and bleeding pump

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



CHARACTERISTICS:

GENERAL:

- ☒ Adult male torso
- ☒ Skin safe prosthetic grade silicone
- ☒ Includes bleeding tubes, bleeding pump
- ☒ Robust, lightweight, washable & reusable
- ☒ Dimensions: 60 × 28 × 79 cm
- ☒ Weight: 14kg

FEATURES:

- ☒ This packable wound trainer simulates the procedures of intraosseous insertions, needle decompression and wound packing.
- ☒ The realistic wound is located to either the left or right of the pelvis, and can be filled with blood via an internal bleeding tube.
- ☒ The two wounds include a rear gunshot entry, and a deep frontal exit wound with a large internal cavity. These features also include built-in tubes with bleeding capabilities.
- ☒ Due to the handmade nature of SIMBODIES products, each trainer can vary slightly from the associated image.
- ☒ We recommend using our Floor Blood Concentrate mixed with water to simulate blood.

INCLUDE:

- ☒ 1 Trauma Torso Task Trainer
- ☒ 1 Manual Blood pump
- ☒ 1 Simulated blood

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>