

HAL AUSCULTATION, AIRWAY & CPR SKILLS TRAINER - OMNI2

PRODUCT NUMBER: 1220



DESCRIPTION:

HAL® is an easy-to-use portable airway management and CPR skills trainer designed for teaching resuscitation and auscultation skills to students and professionals in EMS, Medicine, and Nursing. HAL® features a high-fidelity airway, realistic heart and lungs sounds, and real-time CPR quality feedback to help you teach effective cardiopulmonary resuscitation and airway management.

HAL® is also available in a full-body configuration.

SKILLS:

☑- Airway management

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- ☒- Nasal and oral intubation
- ☒- Resuscitation training (CPR)
- ☒- BVM ventilation
- ☒- Cardiac and pulmonary auscultation assessment
- ☒- Debriefing tool

CHARACTERISTICS:

FEATURES:

- ☒- Adult-sized, male upper torso
- ☒- Anatomically correct landmarks for proper CPR hand placement

AIRWAY MANAGEMENT:

- ☒- Realistic airway with tongue, vocal cords, and trachea
- ☒- Realistic head tilt/chin lift and jaw thrust
- ☒- Nasal and oral Intubation with standard adjuncts: ETT, LMA, King LT®
- ☒- CPR hand position landmarks
- ☒- Visible chest rise with bag-valve-mask ventilation
- ☒- Esophageal intubation and/or excessive ventilation produce gastric distension
- ☒- Right mainstem intubation presents unilateral chest rise

CPR:

- ☒- Compliant with 2015 adult CPR guidelines
- ☒- Realistic chest cavity resistance and recoil
- ☒- Real-time feedback of CPR quality: Compression depth and rate, Ventilation Rate, Excessive Ventilation, No-Flow Time and CPR Cycles

CARDIAC/CIRCULATION:

- ☒- Palpable carotid pulses – Manually generated
- ☒- 40 Programmable abnormal and normal heart and lung sounds

INCLUDE:

- ☒- 1 HAL Skills trainer
- ☒- 1 OMNI2 Wireless Control Interface (Tablet+License)
- ☒- 1 Carrying case
- ☒- 1 Power supply adapter
- ☒- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>