

# SIMBODIES - Adult Standard Female simulator

PRODUCT NUMBER: 1333



## DESCRIPTION:

The **SIMBODIES Advanced Female Manikins** are life-like patient simulators, often used by educators and clinicians for medical training purposes. Their look and feel resembles that of genuine human skin, colouring, hair and eyes, all of which enable the student or trainee to experience a truly immersive training experience.

They are produced using skin safe Prosthetic Grade Platinum Silicone, and are able to **feature medical capabilities** using the replication of **veins, airways, jaws and teeth for added accuracy.**

## MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

## CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>

## SKILLS:

- ☒ Naso and oropharyngeal intubation
- ☒ Supraglottic airway device insertion
- ☒ Laryngoscopy and tracheal intubation
- ☒ Suction skills training
- ☒ Nostrils for nasogastric tube placement
- ☒ Postural airway management
- ☒ BLS and CPR training
- ☒ Application of AED or Defibrillator pads
- ☒ Artificial ventilation exercises: direct, with BVM, ventilator, etc.
- ☒ Patient mobilization and immobilization
- ☒ Direct and indirect pressure application
- ☒ Tourniquet application
- ☒ Pelvic binder, traction device, manual inline stabilisation and spinal collar application
- ☒ Wound packing exercises
- ☒ Circulation assessment
- ☒ Airway sounds auscultation
- ☒ Intraosseous insertions in bilateral humeral head & tibia
- ☒ Needle decompression in the 2nd intercostal space in the midclavicular line (front)
- ☒ Application of chest seal
- ☒ Patient care: bed bathing, undressing, temperature take, etc.

## CHARACTERISTICS:

### GENERAL:

- ☒ Human hair wig
- ☒ Splashproof and washable skin
- ☒ Weight: 25-30kg
- ☒ Age: 20-30 years old
- ☒ Clothing size: M
- ☒ Height: 168cm
- ☒ Ethnic skin tones available

### PATIENT CARE:

- ☒ Head, chin and jaw thrust
- ☒ Mouth and jaw capable or airway management
- ☒ Able to vomit or make vomit noise
- ☒ Close eyelids capable of manual opening
- ☒ Realistic male genitalia

### MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



- ☒ Integrated bluetooth speaker to simulate voice and breathing sounds
- ☒ Fully flexible steel frame (including front reinforced flexible rubber ribcage) for Basic Life Support exercise

#### **INJECTION AND INFUSION:**

- ☒ Various intraosseous injection sites
- ☒ Add IV and venous therapy capabilities with its easy-to-use SIMSLEEVES

#### **TRAUMA:**

- ☒ Needle decompression into the 2nd intercostal space in the midclavicular line
- ☒ Able to bleed at multiple/different sites
- ☒ Emergency devices can be used with the manikin: tourniquet, Intraosseous injection, hemostatic pads, etc.
- ☒ Realistic limbs for patient immobilization techniques: spinal collar, manual inline stabilisation, pelvic binder, etc.
- ☒ Different intraosseous access points
- ☒ Wound packing capabilities (option)

### **INCLUDE:**

1 SIMBODIES adult manikin

### **ADDITIONAL MODULES:**

- ☒ Female genitalia for catheterization (option)
- ☒ Total waterproof capabilities (option)
- ☒ Detachable amputated arm and/or foot (option)
- ☒ Pregnancy 7-8 months suit (option)
- ☒ Needle thoracocentesis (option)
- ☒ SIMSLEEVES - Cannulation arm (option)
- ☒ Rugged transport case (option)

#### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

#### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>