

ZACK - Multipurpose male care

PRODUCT NUMBER: 1198



DESCRIPTION:

The ZACK™ Multipurpose Male Care Simulator is a **full-sized, adult male lower torso** (abdomen and pelvis) designed to provide health care professionals with hands-on **training in the skills required to perform men's healthcare procedures**.

It is useful for demonstrating and practicing catheterization exercises, prostate examinations, testicular examinations, no-scalpel vasectomy, and rectal/colon examinations.

SKILLS:

- ☑- Catheterisation exercises
- ☑- Prostate examination
- ☑- Testicular examination

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- ☒- No scalpel vasectomy
- ☒- Rectum and colon examinations

CASES AND PATHOLOGIES:

- ☒- Melanoma
- ☒- Ulcerative adenocarcinoma
- ☒- Adenomatous polyp
- ☒- Sessile adenocarcinoma
- ☒- Advanced polypoid adenoma
- ☒- Pedunculated adenomatous polyp
- ☒- Pedunculated polyp

CHARACTERISTICS:

- Full-size male lower torso
- Internal bladder for catheterization
- Four interchangeable prostates depicting a moderately enlarged benign prostate, a prostate with two discreet nodules, a prostate with an easily palpable large mass, and a prostate with malignant cancer
- Realistic penis and two scrotal sacks. One scrotal sack is normal, and the other contains tumors in each testicle
- No-Scalpel Vasectomy kit containing two removable scrotal skins, two testicles, and two long vas assemblies that can be advanced as needed for NSV exercises
- Rectum and colon containing benign and malignant masses easily visualized using an appropriate endoscope

INCLUDE:

- ☒- 1 Lower torso male simulator
- ☒- 1 Internal bladder for catheterization
- ☒- 4 Interchangeable prostates
- ☒- 1 No scalpel vasectomy kit
- ☒- 1 Rectum and colon with benign and malignant masses
- ☒- Talcum powder
- ☒- 1 Soft carrying bag
- ☒- 1 Instructions manual

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>