

SIMBODIES - Child standard male simulator

PRODUCT NUMBER: 1334



DESCRIPTION:

The function of **SIMBODIES Standard Child Manikin** is to replicate the look and appearance of a **male infant patient**, around the age of 10 years old, and are produced using skin safe Prosthetic Grade Platinum Silicone.

They have the **appearance and feel of genuine child human skin**, pigmentation, hair and eyes, and are available in both male and female, helping trainees to experience a highly realistic and accurate training experience for all users.

SKILLS:

- ☑- Naso and oropharyngeal intubation
- ☑- Supraglottic airway device insertion
- ☑- Laryngoscopy and tracheal intubation

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- ☒- Suction skills training
- ☒- Nostrils for nasogastric tube placement
- ☒- Postural airway management
- ☒- BLS and CPR training
- ☒- Application of AED or Defibrillator pads
- ☒- Artificial ventilation exercises: direct, with BVM, ventilator, etc.
- ☒- Patient mobilization and immobilization
- ☒- Direct and indirect pressure application
- ☒- Tourniquet application
- ☒- Pelvic binder, traction device, manual inline stabilisation and spinal collar application
- ☒- Wound packing exercises
- ☒- Circulation assessment
- ☒- Airway sounds auscultation
- ☒- Intraosseous insertions in bilateral humeral head & tibia
- ☒- Needle decompression in the 2nd intercostal space in the midclavicular line (front)
- ☒- Application of chest seal
- ☒- Patient care: bed bathing, undressing, temperature take, etc.

CHARACTERISTICS:

GENERAL:

- ☒- Human hair wig
- ☒- Splashproof and washable skin
- ☒- Weight: 15kg
- ☒- Age: 9-10 years old
- ☒- Clothing size: aged 9/10
- ☒- Height: 128cm
- ☒- Ethnic skin tones available

PATIENT CARE:

- ☒- Head, chin and jaw thrust
- ☒- Mouth and jaw capable of airway management
- ☒- Able to vomit or make vomit noise
- ☒- Close eyelids capable of manual opening
- ☒- Realistic male genitalia
- ☒- Integrated bluetooth speaker to simulate voice and breathing sounds
- ☒- Fully flexible steel frame (including front reinforced flexible rubber ribcage) for Basic Life Support exercise

INJECTION AND INFUSION:

- ☒- Various intraosseous injection sites
- ☒- Add IV and venous therapy capabilities with its easy-to-use SIMSLEEVES

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



TRAUMA:

- ☒- Needle decompression into the 2nd intercostal space in the midclavicular line
- ☒- Able to bleed at multiple/different sites
- ☒- Emergency devices can be used with the manikin: tourniquet, Intraosseous injection, hemostatic pads, etc.
- ☒- Realistic limbs for patient immobilization techniques: spinal collar, manual inline stabilisation, pelvic binder, etc.
- ☒- Different intraosseous access points
- ☒- Wound packing capabilities (option)

INCLUDE:

- ☒- 1 SIMBODIES child manikin
- ☒- 1 IO leg pod (installed)
- ☒- 1 IO humeral pod (installed)
- ☒- 1 Realistic and intubatable airway (installed)

ADDITIONAL MODULES:

- ☒- Total waterproof capability (option)
- ☒- Fixed amputated arm or foot (option)
- ☒- Detachable amputated arm or foot (option)
- ☒- Integrated bleeding tubes for hemorrhage control (option)
- ☒- Catheterization capability (option)
- ☒- Rugged transport case (option)
- ☒- Blood glucose testing (option)
- ☒- Fractured pelvis (option)
- ☒- Surgical airway capability (option)
- ☒- SIMSLEEVES – Cannulation sleeve arm (option)

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>