



TRUPICC - Ultrasound guided IV and PICC simulator

PRODUCT NUMBER: 1130



DESCRIPTION:

TruPICC is an adult male training arm designed for effective training in **peripherally inserted central catheters and IV insertion**.

This **IV and PICC training model** develops the user's skills associated with needle placement, guidewires and catheters using ultrasound. It's also ideal for ultrasound manufacturers for education and demonstrations.

SKILLS:

- ☑- Probe positioning and movement
- ☑- Recognition of arterial and vein vessels and nerves in soft responsive tissue
- ☑- Using ultrasound to target a nerve for ultrasound guided regional anesthesia
- ☑- Ability to verify needle tip location and to practice the entire regional anesthesia procedure

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



- ☒- Visualizing of the artery and vein laterally beside the nerve
- ☒- IV training
- ☒- PICC line training
- ☒- Fracture assessment detection of a fracture and diagnosing bone stress injury

CHARACTERISTICS:

- ☒- The upper arm anatomy features the brachial vein, basilic vein, cephalic, median cubital vein, superior vena cava and the brachial artery for vessel differentiation.
- ☒- Differentiation between the basilic vein and the basilic artery.
- ☒- Realistic blood flashback upon entry into the vessels. Verify correct needle placement by withdrawal of fluids.
- ☒- Option to switch on/off flow to either the cephalic or basilic veins to prolong the life of the insert.
- ☒- Real feel vascular 'tenting' upon entry into the vessel
- ☒- Dimensions: 90cm x 40cm x 20cm
- ☒- Weight: 6kg aprox.

INCLUDE:

- ☒- 1 TruPICC model ready to use
- ☒- 1 TruPICC carrier case
- ☒- 1 USB pen drive user manual
- ☒- 1 Bottle of artificial blood concentrate (250mL)
- ☒- 1 Luerlock syringe for fluidsystem preparation/withdrawal

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>