

TruBaby X - Multipurpose Pediatric Simulator

PRODUCT NUMBER: 1131



DESCRIPTION:

Our **pediatric clinical skills training** model is ideal for anesthetists, nurses & other pediatric emergency medical professionals practicing Directly Observed Practical Skills (DOPS), Pediatric Advanced Life Support (PALS) and critical emergency medicine including infant **airway management, CPR, IV cannulation, peripheral venous cannulation, PICC line insertion, needle thoracentesis & chest drain, urethral catheterization and lumbar puncture.**

TruBaby X is incredibly lifelike with the appearance, weight, size and movement of a 5-month-old infant.

Can be used along with our patient monitor simulator (option) to integrate clinical decision-making with medical technique training.

SKILLS:

☑- Airway management

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>

- ☒- Peripheral venous cannulation (arm, hand and foot)
- ☒- Lumbar puncture
- ☒- PICC line insertion
- ☒- Ultrasound guided PICC line insertion
- ☒- Needle thoracentesis for a tension pneumothorax
- ☒- Chest drain (air only)
- ☒- IO infusion tibia
- ☒- Urethral catheterization
- ☒- CPR

CHARACTERISTICS:

- ☒- Realistic & durable AirSim X airway with 5-year warranty
- ☒- Realistic movement including head tilt, chin lift & jaw thrust
- ☒- Interchangeable head in case of repairs
- ☒- Realistic rib structures including xiphoid process & clavicle
- ☒- Palpable landmarks in the 2nd intercostal space midclavicular line & 5th intercostal space mid-axillary line
- ☒- Palpable vertebrae landmarks
- ☒- Lumbar puncture can be practiced at L3-L4 & L4-L5 spine locations
- ☒- Proximal tibia interosseous needle insertion site
- ☒- Contains tibial tuberosity anatomy & patella anatomy
- ☒- Realistic anatomy to practice inserting catheter tube to extract urine
- ☒- Option for interchangeable male and female genitalia
- ☒- Optional arm facilitates ultrasound-guided IV & PICC line insertion
- ☒- Fully enclosed fluid management system giving realistic blood flashback & flow
- ☒- Allows needle cannulation at various various locations in hand, arm & foot with realistic blood flashback
- ☒- Realistic look & feel of tissues
- ☒- Lifelike responsiveness during procedure training
- ☒- Set-up time is less than 5 minutes

INCLUDE:

- ☒- 1 TruBaby X model
- ☒- 1 Syringewith yellow and white tubing and quick fit connectors
- ☒- 4 Needle Decompression inserts (2 Left & 2 Right)
- ☒- 4 Chest Drain inserts (2 Left & 2 Right)
- ☒- 1 TruBaby X Ultrasound Arm (Left)
- ☒- 1 TruBaby X Ultrasound Arm (Right)
- ☒- 10 TruBaby X IO inserts (5 Left & 5 Right)
- ☒- 1 TruBaby X Lumbar Puncture insert
- ☒- 1 Trubaby X Male Genitalia insert

MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>



More Than Simulators™
Excellence. Every patient. Every time.

- ☒- 1 Bottle of artificial blood concentrate(250ml)
- ☒- 1 BottleTruCorp Lubrication (100mL)
- ☒- 1 USB pen drive user manual
- ☒- 1 TruBaby X carrier case



MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º
44570 Calanda (Teruel)
España

CONTACTO

sales@morethansimulators.com
(+34) 936 06 77 05 / (+351) 21 123 85 18
<https://morethansimulators.com>