



# Brayden Basic

PRODUCT NUMBER: 1657



## DESCRIPTION:

**Life-like adult torso for CPR training**, fully compliant with AHA and ERC (2015) Guidelines

## SKILLS:

### CHEST COMPRESSION:

☑- Chest resistance is representative of an average adult. Audible feedback through 'clicker' reinforces correct compression depth (5cm).

### VENTILATION:

☑- Ventilation can be achieved through Mouth to Mouth, pocket mas or Bag-Valve-Mask.

## MORE THAN SIMULATORS S.L.

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

## CONTACTO

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>



## CHARACTERISTICS:

- ☒ Oral and nasal passages which allow realistic nose pinch.
- ☒ Highly realistic movement with correct ventilations.
- ☒ Head tilt and chin lift for opening airway.
- ☒ Anatomically correct, harmless and latex-free skin.
- ☒ Easy AED training pad removal.
- ☒ RoHS and REACH compliant.
- ☒ Air vent reduces the risk of cross contamination. Air from Brayden's disposable lungs exits through the back of his head.

## INCLUDE:

- ☒ 1 Brayden Manikin
- ☒ 1 Mat Bag
- ☒ 1 Extra Face Skin
- ☒ 10 Face Shields
- ☒ 1 Artificial Lung
- ☒ 1 Instructions for Use

### **MORE THAN SIMULATORS S.L.**

Calle de la Hoya, 4 2º  
44570 Calanda (Teruel)  
España

### **CONTACTO**

sales@morethansimulators.com  
(+34) 936 06 77 05 / (+351) 21 123 85 18  
<https://morethansimulators.com>